


# RNA 6000 Ladder

Store below  $-70^{\circ}\text{C}$ .

<b>Catalog # (P/N):</b>	AM7152
<b>Volume:</b>	3 x 20 $\mu\text{L}$
<b>Concentration:</b>	150 ng/ $\mu\text{L}$
<b>Product Description:</b>	The RNA 6000 Ladder is a set of six transcripts, 0.2, 0.5, 1.0, 2.0, 4.0 and 6.0 kb in length, designed for use with the Agilent 2100 Bioanalyzer. The 0.2 kb RNA is at a concentration of 20 ng/ $\mu\text{L}$ ( $A_{260}$ ).
<b>Storage Conditions:</b>	Store below $-70^{\circ}\text{C}$ . Avoid multiple freeze-thaw cycles. The product may be stored short-term at $-20^{\circ}\text{C}$ .
<b>Storage Buffer:</b>	0.1 mM EDTA

## USER INFORMATION

**General Information:** Slight variations in the signal intensities of the ladder peaks may be seen between different reagent lots (as shown in Figure 1). These minor variations will not compromise assay performance.



**Figure 1:** Typical signal intensity variations seen between two different lots of RNA 6000 Ladder.

**Handling Instructions:**

This product is very sensitive to degradation by exogenous ribonucleases introduced during handling. Wear gloves when handling this product. Use RNase-free reagents, tubes, and barrier pipette tips. **A degraded ladder can lead to erroneous quantitation.**

**Thawing instructions**

Thaw just to completion at 37°C, vortex for a few seconds when fully thawed, and place on ice. Aliquot the RNA, if necessary, to minimize freeze-thaw cycles ( $\leq 5$ ).

**Applications:**

The RNA 6000 Ladder is designed for use with the Agilent 2100 bioanalyzer and the RNA 6000 LabChip® kits.

Prior to use, centrifuge the tube for a few seconds to collect the contents of the bottom of the tube. Each of the three tubes contains 20  $\mu\text{L}$  of RNA 6000 Ladder. Since only 1  $\mu\text{L}$  is used per chip, it may be advisable to aliquot the product into several RNase-free microcentrifuge tubes. Generally when degradation occurs, the high molecular weight peaks disappear. If the RNA peak pattern changes dramatically from the examples on this data sheet, the product may be degraded and a fresh tube of RNA 6000 Ladder should be used.

**Ladder Denaturation**

The RNA 6000 ladder must be heat denatured before use. Each day, dispense a 1-day supply of the product into a fresh microcentrifuge tube and heat the tube at 70°C for 2 min. Place the tube directly on ice for 5 min to snap-cool. Briefly centrifuge to clear any condensate from the tube's walls and cap. Keep the RNA 6000 Ladder on ice throughout the day.

---

**QUALITY CONTROL**

---

**Functional Testing:**

150 ng (1  $\mu\text{L}$ ) will generate 6 distinct bands when run on a RNA LabChip using the Total RNA or mRNA assay. The ladder remains intact when incubated at 37°C overnight.

---

**OTHER INFORMATION**

---

**Material Safety Data Sheets:**

Material Safety Data Sheets (MSDSs) can be printed or downloaded from product-specific links on our website at the following address: [www.ambion.com/techlib/msds](http://www.ambion.com/techlib/msds). Alternatively, e-mail your request to [MSDS\\_Inquiry\\_CCRM@appliedbiosystems.com](mailto:MSDS_Inquiry_CCRM@appliedbiosystems.com). Specify the catalog or part number(s) of the product(s), and we will e-mail the associated MSDSs unless you specify a preference for fax delivery. For customers without access to the internet or fax, our technical service department can fulfill MSDS requests placed by telephone or postal mail. (Requests for postal delivery require 1–2 weeks for processing.)

**Warranty and Liability:**

***For research use only. Not for use in diagnostic procedures.***

Applied Biosystems is committed to delivering superior product quality and performance, supported by industry-leading global service and technical support teams. Warranty information for the accompanying consumable product is available at [www.ambion.com/info/warranty](http://www.ambion.com/info/warranty) in "Limited Warranty for Consumables," which is subject to the exclusions, conditions, exceptions, and limitations set forth under the caption "EXCLUSIONS, CONDITIONS, EXCEPTIONS, AND LIMITATIONS" in the full warranty statement. Please contact Applied Biosystems if you have any questions about our warranties or would like information about post-warranty support.

Information in this document is subject to change without notice. Applied Biosystems assumes no responsibility for any errors that may appear in this document.

Applied Biosystems disclaims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall Applied Biosystems be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

**Trademarks, Patents, and Licensing:**

Applied Biosystems, AB (Design), and Ambion are registered trademarks of Applera Corporation or its subsidiaries in the US and/or certain other countries. LabChip is a US trademark of Caliper Technologies Corp. All other trademarks are the sole property of their respective owners.

© 2008 Ambion, Inc. All rights reserved. 4382133B